

# AMERCARE

PROTECTING YOUR PRODUCT AND ENVIRONMENT

## Easy Change Cuff Ring System



Amercure's Easy Change Cuff Ring system has been specifically designed to allow a sterile glove change to be carried providing operator protection.

Whether you need to carry out a routine change, replace a damaged glove simply change the glove size the Easy Change system enables you to change gloves quickly and easily without breaching the integrity of your environment.

The Cuff Ring utilises a double lip seal arrangement that ensures that the containment is never broken.

The glove is changed by loading a sterile glove onto an inner ring inside the Isolator and then simply pushing the other used glove and ring through into the sleeve or into a bag.



## Easy Change Cuff Ring

- Double lip seal ensures containment is maintained throughout change procedure.
- All parts of system are suitable for sterilisation including autoclaving at 121 C.
- Can be used with most makes of beaded glove.
- Allows an operator safe change as well as a sterile change.
- Provision for bagging out of contaminated gloves.
- Compatible with a wide range of cleaning and sterilising agents.
- Can be used with any make of Isolator including flexible film Isolators and Half Suits.
- 90mm internal diameter gives comfortable clearance over hand and wrist

AMERCARE LTD  
UNIT L, OYO BUSINESS UNITS  
HOWLAND ROAD  
THAME, OXFORDSHIRE UK  
OX9 3GQ

TELEPHONE  
+44 (0)1844 261244  
E-MAIL  
sales@amercare.co.uk

FACSIMILE  
+44 (0)1844 261611

WEB SITE  
www.amercare.co.uk

