

AMERCARE

PROTECTING YOUR PRODUCT AND ENVIRONMENT

ELS Neoprene Gloves



Amercure is a market leader in the design, manufacture and supply of rigid isolators for Pharmacy and Nuclear Medicine Applications and their consumable spares.

All our sleeves, gauntlets and gloves give a high level of control and comfort to give you the best manipulation of your product.

All the materials are to the highest quality and made to our exacting standards.

Neoprene ELS Gloves are strong and hard wearing.

The material allows for a good tight seal onto the either the Amercare Easy Change Cuff System or our Oval Cuffs.

The glove is very resistant to both chemical and physical wear.

We offer this style of Glove in Hand Sizes 7 and 8 from stock
Size 9's are available as special order items

Amercure Code	Description	Size
ISN 4360	ELS Neoprene Glove	Size 7
ISN 4361	ELS Neoprene Glove	Size 8